

User Fitting Form

www.yourkoalaa.com

Step 1. Take Personal Details

Users Full Name						
Email			Phone Number			
Fitting Location			Date Of Birth			
Gender Male	Female		Pre	Prefer not to say		\bigcirc
A A	B					
Above Elbow	Below Elbow	\bigcirc	Through Wrist	\bigcirc	Partial Hand	\bigcirc
Limb History	Amputation	\bigcirc	Congenital	\bigcirc		
Affected Side	Left	\bigcirc	Right	\bigcirc		
Clinicians Details						
Clinicians Name			Clinic Name			
Email			Phone Numb	er		

Step 2. What would you like to do with your Koalaa?

Write down the main activity or task the user would like to achieve with their Koalaa prosthetic

How do you currently find this activity?


Circle the answer that relates to the users current experience with the activity





Step 3. Take Measurements


Use the below guide for below elbow or through wrist limb difference




Measure the circumference of the widest part of limb

Α

С



Measure the circumference of the narrowest part of limb



Measure length from crook of elbow to end of limb



cm

В

Ν

Measure length of opposite arm from crook of elbow to middle of palm

cm

No

cm

Use the below guide for partial hand or partial finger wrist limb difference



quotes and stories) publicly and with our partners?



Step 4. Fit the Koalaa and document items given to the user

Sleeve Type

